



CAT. NO. SS	CAT. NO. TT	1/3RD TUBULAR PLATES FOR 3.5 MM SCREWS
291.004	291.044	1/3RD TUBULAR PLATES FOR 3.5 MM SCREWS - 04 HOLES, L - 52MM
291.005	291.055	1/3RD TUBULAR PLATES FOR 3.5 MM SCREWS - 05 HOLES, L - 64MM
291.006	291.066	1/3RD TUBULAR PLATES FOR 3.5 MM SCREWS - 06 HOLES, L - 76MM
291.007	291.077	1/3RD TUBULAR PLATES FOR 3.5 MM SCREWS - 07 HOLES, L - 88MM
291.008	291.088	1/3RD TUBULAR PLATES FOR 3.5 MM SCREWS - 08 HOLES, L - 100MM
291.009	291.099	1/3RD TUBULAR PLATES FOR 3.5 MM SCREWS - 09 HOLES, L - 112MM
291.010	291.100	1/3RD TUBULAR PLATES FOR 3.5 MM SCREWS - 10 HOLES, L - 124MM
291.011	291.111	1/3RD TUBULAR PLATES FOR 3.5 MM SCREWS - 11 HOLES, L - 136MM
291.012	291.120	1/3RD TUBULAR PLATES FOR 3.5 MM SCREWS - 12 HOLES, L - 148MM

Specification

- One-Third Tubular plate is available in both Stainless Steel and Titanium material.
- The plate holes support 3.5mm cortex and 4.0mm cancellous screws.

Indication for use

- One Third Tubular Plate is used to fix Lateral Malleolar fracture, Distal Ulnar fracture and Olecranon fracture. The plate neutralizes any shearing or bending force which could disrupt an unprotected lag screw fixation.

