


3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION



CAT. NO. SS	CAT. NO. TT	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION
348.523L	648.523L	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 3H LEFT
348.523R	648.523R	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 3H RIGHT
348.524L	648.524L	3.5MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 4 HOLE LEFT
348.524R	648.524R	3.5MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 4 HOLE RIGHT
348.525L	648.525L	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 5H LEFT
348.525R	648.525R	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 5H RIGHT
348.526L	648.526L	3.5MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 6 HOLE LEFT
348.526R	648.526R	3.5MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 6 HOLE RIGHT
348.527L	648.527L	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 7H LEFT
348.527R	648.527R	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 7H RIGHT
348.528L	648.528L	3.5MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 8 HOLE LEFT
348.528R	648.528R	3.5MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 8 HOLE RIGHT
348.529L	648.529L	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 9H LEFT
348.529R	648.529R	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 9H RIGHT
348.530L	-	3.5MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 10 HOLE LEFT
348.530R	-	3.5MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 10 HOLE RIGHT
348.531L	-	3.5MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 11 HOLE LEFT
348.531R	-	3.5MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITHOUT LATERAL EXTENSION 11 HOLE RIGHT

Specification

- Plates (Right & Left) available in both – Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 2.7-3.5mm cortical and Wise Lock Screws.
- Plate available with 3,4,5,6,7,8,9,10,11 holes.

Indication for use

These plates are indicated for:

- Intraarticular fractures of the distal humerus
- Supracondylar fractures of the distal humerus
- Non-unions of the distal humerus

