


3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITH LATERAL EXTENSION



CAT. NO. SS	CAT. NO. TT	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITH LATERAL EXTENSION
348.453L	648.453L	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITH LATERAL EXTENSION 3H LEFT L 78
348.453R	648.453R	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITH LATERAL EXTENSION 3H RIGHT L 78
348.455L	648.455L	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITH LATERAL EXTENSION 5H LEFT L 104
348.455R	648.455R	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITH LATERAL EXTENSION 5H RIGHT L 104
348.457L	648.457L	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITH LATERAL EXTENSION 7H LEFT L 130
348.457R	648.457R	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITH LATERAL EXTENSION 7H RIGHT L 130
348.459L	648.459L	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITH LATERAL EXTENSION 9H LEFT L 156
348.459R	648.459R	3.5 MM LCP DORSOLATERAL LOWER HUMERUS PLATE WITH LATERAL EXTENSION 9H RIGHT L 156

Specification

- Plates (Right & Left) available in both – Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 2.7-3.5mm cortical and Wise Lock Screws.
- Plate available with 3,5,7 & 9 holes

Indication for use

- Intraarticular fractures of the distal humerus
- Supracondylar fractures of the distal humerus
- Non-unions of the distal humerus