

## 4.5MM CURVED BROAD LCP PLATES



CAT. NO. SS	CAT. NO. TT	4.5MM Curved Broad LCP Plates
310.230	610.230	4.5MM Curved Broad LCP Plates 10 holes
310.232	610.232	4.5MM Curved Broad LCP Plates 12 holes
310.234	610.234	4.5MM Curved Broad LCP Plates 14 holes
310.236	610.236	4.5MM Curved Broad LCP Plates 16 holes
310.238	610.238	4.5MM Curved Broad LCP Plates 18 holes
310.240	610.240	4.5MM Curved Broad LCP Plates 20 holes

### Specification

- Available in both Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5mm/5.0mm cortical and Wise-lock screw.
- Plate available with 10,12,14,16,18 & 20 holes.

### Indication for use

- It is used to fix the Femoral shaft fracture.
- The curve in the plate easily fit over the Femur.

