



CAT. NO. SS	CAT. NO. TT	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH.
229.135	629.135	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 0 35 MM
229.140	629.140	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 040 MM
229.150	629.150	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 050 MM
229.152	629.152	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 052 MM
229.154	629.154	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 054 MM
229.155	629.155	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 055 MM
229.156	629.156	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 056 MM
229.160	629.160	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 060 MM
229.165	629.165	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 065 MM
229.170	629.170	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 070 MM
229.175	629.175	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 075 MM
229.180	629.180	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 080 MM
229.185	629.185	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 085 MM
229.190	629.190	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 090 MM
229.195	629.195	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 095 MM
229.200	629.200	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 100 MM

Specification

- Available in 35mm to 110mm length with 6.5mm diameter
- 16mm threaded shaft
- Available in both Stainless Steel & Titanium

Indication for use

- Cannulated screws are indicated for securing plates into metaphyseal and epiphyseal areas of bones
- Deep cuts and widely spaced threads enhances the holding power of the screw in soft bones