

LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH.



CAT. NO. SS	CAT. NO. TT	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH.
229.135	629.135	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 0 35MM
229.140	629.136	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 040 MM
229.145	629.137	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 045 MM
229.150	629.138	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 050 MM
229.152		LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 052 MM
229.154		LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 054 MM
229.155	629.141	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 055 MM
229.156		LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 056 MM
229.160	629.143	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 060 MM
229.165	629.144	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 065 MM
229.170	629.145	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 070 MM
229.175	629.146	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 075 MM
229.180	629.147	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 080 MM
229.185	629.148	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 085 MM
229.190	629.149	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 090 MM
229.195	629.150	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 095 MM
229.200	629.151	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 100 MM
229.205	629.152	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 105 MM
229.210	629.153	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH. X 1 10 MM

Specification

- Available in 35mm to 110mm length with 6.5mm diameter
- Partially 16mm threaded shaft
- Available in both Stainless Steel & Titanium

Indication for use

- 6.5mm Cannulated screws are indicated for securing plates into metaphyseal and epiphyseal areas of bones
- Deep cuts and widely spaced threads enhances the holding power of the screw in soft bones