



CAT. NO. SS	CAT.NO. TT	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH.
229.235	649.172	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 0 35 MM
229.240	649.174	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 040 MM
229.245	649.176	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 045 MM
229.250	649.178	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 050 MM
229.252		LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 052 MM
229.254		LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 054 MM
229.255	649.184	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 055 MM
229.256		LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 056 MM
229.258		LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 058 MM
229.260	649.188	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 060 MM
229.265	649.190	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 065 MM
229.270	649.192	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 070 MM
229.275	649.194	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 075 MM
229.280	649.196	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 080 MM
229.285	649.198	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 085 MM
229.290	649.200	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 090 MM
229.295	649.202	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 095 MM
229.300	649.204	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 100 MM
229.305	649.206	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 105 MM
229.310	649.208	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH. X 1 10 MM

## Specification

- Available in 35mm to 110mm length with 6.5mm diameter
- Partially 32mm threaded shaft
- Available in both Stainless Steel & Titanium

## Indication for use

- Cannulated screws are indicated for securing plates into metaphyseal and epiphyseal areas of bones
- Deep cuts and widely spaced threads enhances the holding power of the screw in soft bones