



CAT. NO. SS	CAT. NO. TT	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.)
229.025	629.425	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.) X 0 25 MM
229.040	629.440	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.) X 040 MM
229.045	629.445	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.) X 045 MM
229.050	629.450	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.) X 050 MM
229.055	629.455	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.) X 055 MM
229.060	629.460	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.) X 060 MM
229.065	629.465	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.) X 065 MM
229.070	629.470	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.) X 070 MM
229.075	629.475	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.) X 075 MM

Specification

- Compress the plate to the bone and provide inter fragmentary compression.
- Available in 30 mm to 75 mm length
- Fully threaded shaft.
- Self-tapping tip.
- Available in both Titanium and stainless steel.

Indication for use

- It is used in the fixation of Large Fragment Wise-Lock Plates.

