



CAT. NO. SS	CAT. NO. TT	ANGULAR STABLE LOCKING SCREW (ASLS)
268.828	668.828	Angular Stable Locking Screw (ASLS) 4.9mm x 28mm
268.830	668.830	Angular Stable Locking Screw (ASLS) 4.9mm x 30mm
268.832	668.832	Angular Stable Locking Screw (ASLS) 4.9mm x 32mm
268.834	668.834	Angular Stable Locking Screw (ASLS) 4.9mm x 34mm
268.836	668.836	Angular Stable Locking Screw (ASLS) 4.9mm x 36mm
268.838	668.838	Angular Stable Locking Screw (ASLS) 4.9mm x 38mm
268.840	668.840	Angular Stable Locking Screw (ASLS) 4.9mm x 40mm
268.842	668.842	Angular Stable Locking Screw (ASLS) 4.9mm x 42mm
268.844	668.844	Angular Stable Locking Screw (ASLS) 4.9mm x 44mm
268.846	668.846	Angular Stable Locking Screw (ASLS) 4.9mm x 46mm
268.848	668.848	Angular Stable Locking Screw (ASLS) 4.9mm x 48mm
268.850	668.850	Angular Stable Locking Screw (ASLS) 4.9mm x 50mm
268.855	668.855	Angular Stable Locking Screw (ASLS) 4.9mm x 55mm
268.860	668.860	Angular Stable Locking Screw (ASLS) 4.9mm x 60mm

Specification

- Available in both Stainless steel & Titanium.
- Available in 28mm to 60mm length with 4.9mm diameter.

Indication for use

- Angular Stable Locking Screw Provides improved stiffness and decreased inter fragmentary movement compared to conventional locking.

