



CAT. NO. SS	CAT. NO. TT	BROAD DCP FOR 4.5 MM SCREWS
276.404	676.404	BROAD DCP FOR 4.5 MM SCREWS- 04 HOLES,
276.405	676.405	BROAD DCP FOR 4.5 MM SCREWS- 05 HOLES,
276.406	676.406	BROAD DCP FOR 4.5 MM SCREWS- 06 HOLES, L- 103 MM
276.407	676.407	BROAD DCP FOR 4.5 MM SCREWS- 07 HOLES, L- 119 MM
276.408	676.408	BROAD DCP FOR 4.5 MM SCREWS- 08 HOLES, L- 135 MM
276.409	676.409	BROAD DCP FOR 4.5 MM SCREWS- 09 HOLES, L- 151 MM
276.410	676.410	BROAD DCP FOR 4.5 MM SCREWS- 10 HOLES, L- 167 MM
276.411	676.411	BROAD DCP FOR 4.5 MM SCREWS- 11 HOLES, L- 183 MM
276.412	676.412	BROAD DCP FOR 4.5 MM SCREWS- 12 HOLES, L- 199 MM
276.413	676.413	BROAD DCP FOR 4.5 MM SCREWS- 13 HOLES, L-
276.414	676.414	BROAD DCP FOR 4.5 MM SCREWS- 14 HOLES, L- 231 MM
276.415	676.415	BROAD DCP FOR 4.5 MM SCREWS- 15 HOLES, L-
276.416	676.416	BROAD DCP FOR 4.5 MM SCREWS- 16 HOLES, L- 263 MM
276.417	676.417	BROAD DCP FOR 4.5 MM SCREWS- 17 HOLES, L- 263 MM
276.418	676.418	BROAD DCP FOR 4.5 MM SCREWS- 18 HOLES, L- 263 MM
276.420	676.420	BROAD DCP FOR 4.5 MM SCREWS- 20 HOLES, L- 263 MM

Specification

- Available in both Stainless Steel and Titanium material.
- Plate holes accept 4.5/6.5 mm cortical screws and Cancellous screws
- Plates available with 4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19 & 20 holes.
- Plate bears self compression holes

Indication for use

- Available in both Stainless Steel and Titanium material.
- Plate holes accept 4.5/6.5 mm cortical screws and Cancellous screws
- Plates available with 4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19 & 20 holes.
- Plate bears self compression holes