



CAT. NO. SS	CAT. NO. TT	CALCANEAL PLATE (PLAIN)
282.060	682.060	CALCANEAL PLATE, L - 60MM
282.070	682.070	CALCANEAL PLATE, L - 70MM

Specification

- This plate are available in 60-70mm length.
- All holes of the plate support 3.5mm cortical screws.
- Available in both Titanium & Stainless Steel.

Indication for use

- 3.5mm Calcaneal Plate is used for non-articular body fractures (Non-comminuted and comminuted).