



CAT. NO.	METAL	CALCANEUS ANTERIOR LCP PLATE
618.420L	TT	CALCANEUS ANTERIOR LCP PLATE 5H LEFT
618.420R	TT	CALCANEUS ANTERIOR LCP PLATE 5H RIGHT
618.425R	TT	CALCANEUS POSTERIOR LCP PLATE 5H RIGHT
618.430L	TT	CALCANEUS POSTERIOR LCP PLATE 6H LEFT
618.430R	TT	CALCANEUS POSTERIOR LCP PLATE 6H RIGHT
618.422L	TT	CALCANEUS ANTERIOR LCP LARGE EXTENDED PLATE 5H LEFT
618.422R	TT	CALCANEUS ANTERIOR LCP LARGE EXTENDED PLATE 5H RIGHT

Specification

- Calcaneus Anterior/Posterior LCP Plate is available in Titanium.
- It accepts 3.5 mm Locking Screws with left & right side.

Indication for use

- Locking Calcaneal Plates address complex fractures of the calcaneus.
- The Locking Calcaneal Plate is indicated for fractures and osteotomies of the calcaneus including, but not limited to, extra-articular, intra-articular, joint depression, tongue type, and severely comminuted fractures.