



CAT. NO.	METAL	CALCANEUS COMBINED LCP
618.408L	TT	CALCANEUS COMBINED LCP SMALL EXTENDED 8H Left
618.408R	TT	CALCANEUS COMBINED LCP SMALL EXTENDED 8H Right
618.415L	TT	CALCANEUS COMBINED LCP LARGE EXTENDED 8H Left
618.415R	TT	CALCANEUS COMBINED LCP LARGE EXTENDED 8H Right

Specification

- Calcaneus Combined LCP is available in Titanium.
- It accepts 3.5 mm Locking Screws.

Indication for use

- Locking Calcaneal Plates address complex fractures of the calcaneus.
- The Locking Calcaneal Plate is indicated for fractures and osteotomies of the calcaneus including, but not limited to, extra-articular, intra-articular, joint depression, tongue type, and severely comminuted fractures.