



CAT. NO.	METAL	CANCELLOUS SCREWS Ø 4.5 MM (F.TH)
210.435	SS	CANCELLOUS SCREW 4.5MM F.TH. X 35MM
210.440	SS	CANCELLOUS SCREW 4.5MM F.TH. X 40MM
210.445	SS	CANCELLOUS SCREW 4.5MM F.TH. X 45MM
210.450	SS	CANCELLOUS SCREW 4.5MM F.TH. X 50MM

Specification

- Available in 35,40,45 & 50mm length with 4.5mm diameter.
- Fully threaded shaft.
- Available in Stainless Steel.

Indication for use

- Cancellous Screws are designed specifically for use in the spongy, less dense, and trabecular bone tissue, commonly found at the ends of long bones, such as the femur and the tibia.