



CAT. NO. SS	CAT. NO. TT	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH.
213.025	613.425	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 025 MM
213.030	613.430	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 030 MM
213.035	613.435	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 035 MM
213.040	613.440	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 040 MM
213.045	613.445	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 045 MM
213.050	613.450	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 050 MM
213.055	613.455	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 055 MM
213.060	613.460	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 060 MM
213.065	613.465	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 065 MM
213.070	613.470	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 070 MM
213.075	613.475	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 075 MM
213.080	613.480	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 080 MM
213.085	613.485	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 085 MM
213.090	613.490	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 090 MM
213.095	613.495	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 095 MM
213.100	613.500	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 100 MM
213.105	613.505	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 105 MM
213.110	613.510	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 110 MM
213.115	613.515	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 115 MM
213.120	613.520	CANCELLOUS SCREWS Ø 6.5MM/16 MM TH. X 120 MM

Specification

- Available in 25mm to 120mm length with 6.5mm diameter
- Partially threaded shaft
- Available in both Stainless Steel & Titanium

Indication for use

- Cancellous screws are indicated for securing plates into metaphyseal and epiphyseal areas of bones
- Deep cuts and widely spaced threads enhances the holding power of the screw in soft bones

