



CAT. NO. SS	CAT. NO. TT	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH.
214.030	641.030	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 030MM
214.035	641.035	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 035MM
214.040	641.040	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 040MM
214.045	641.045	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 045 MM
214.050	641.050	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 050 MM
214.055	641.055	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 055 MM
214.060	641.060	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 060 MM
214.065	641.065	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 065 MM
214.070	641.070	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 070 MM
214.075	641.075	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 075 MM
214.080	641.080	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 080 MM
214.085	641.085	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 085 MM
214.090	641.090	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 090 MM
214.095	641.095	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 095 MM
214.100	641.100	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 100 MM
214.105	641.105	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 105 MM
214.110	641.110	CANCELLOUS SCREWS Ø 6.5MM/32 MM TH. X 110 MM

Specification

- Available in 40mm to 120mm length with 6.5mm diameter
- Partially threaded shaft
- Available in both Stainless Steel & Titanium

Indication for use

- Cancellous screws are indicated for securing plates into metaphyseal and epiphyseal areas of bones
- Deep cuts and widely spaced threads enhances the holding power of the screw in soft bones