



CAT. NO. SS	CAT. NO. TT	CANNULATED HERBERT SCREW
230.012	639.012	CANNULATED HERBERT SCREW - , L- 12 MM
230.014	639.014	CANNULATED HERBERT SCREW - , L- 14 MM
230.016	639.016	CANNULATED HERBERT SCREW - , L- 16 MM
230.018	639.018	CANNULATED HERBERT SCREW - , L- 18 MM
230.020	639.020	CANNULATED HERBERT SCREW - , L- 20 MM
230.022	639.022	CANNULATED HERBERT SCREW - , L- 22 MM
230.024	639.024	CANNULATED HERBERT SCREW - , L- 24 MM
230.026	639.026	CANNULATED HERBERT SCREW - , L- 26 MM
230.028	639.028	CANNULATED HERBERT SCREW - , L- 28 MM
230.030	639.030	CANNULATED HERBERT SCREW - , L- 30 MM
	639.032	CANNULATED HERBERT SCREW - , L- 32 MM
	639.034	CANNULATED HERBERT SCREW - , L- 34 MM
	639.035	CANNULATED HERBERT SCREW - , L- 35 MM
	639.036	CANNULATED HERBERT SCREW - , L- 36 MM
	639.038	CANNULATED HERBERT SCREW - , L- 38 MM
	639.040	CANNULATED HERBERT SCREW - , L- 40 MM
	639.042	CANNULATED HERBERT SCREW - , L- 42 MM
	639.044	CANNULATED HERBERT SCREW - , L- 44 MM
	639.045	CANNULATED HERBERT SCREW - , L- 45 MM
	639.046	CANNULATED HERBERT SCREW - , L- 46 MM
	639.048	CANNULATED HERBERT SCREW - , L- 48 MM
	639.050	CANNULATED HERBERT SCREW - , L- 50 MM
	639.055	CANNULATED HERBERT SCREW - , L- 55 MM
	639.060	CANNULATED HERBERT SCREW - , L- 60 MM
	639.065	CANNULATED HERBERT SCREW - , L- 65 MM
	639.070	CANNULATED HERBERT SCREW - , L- 70 MM

Specification

- Herbert Screws have headless design, it means that the screw is completely embedded in the bone, without any protrusions to cause tissue irritation even in intra-articular placement.
- Cannulation helps ensure precise placement of the screw.
- Screws available in both Titanium and Stainless Steel.

Indication for use

- Herbert Screws are designed for the management of fractures where minimal tissue coverage makes standard screw use inappropriate or where extreme precision in fragment alignment is imperative.
- Such locations might include fractures of the: Olecranon and malleolar region

Distal femur

Proximal tibia

Distal tibia