

CANNULATED SCREWS Ø 4.0MM (FULLY THREADED)



CAT. NO. SS	METAL	CANNULATED SCREWS Ø 4.0MM (fully Threaded)
217.210	SS	Cannulated Screws Ø 4.0mm (F.th) x 10mm
217.212	SS	Cannulated Screws Ø 4.0mm (F.th) x 12mm
217.214	SS	Cannulated Screws Ø 4.0mm (F.th) x 14mm
217.216	SS	Cannulated Screws Ø 4.0mm (F.th) x 16mm
217.218	SS	Cannulated Screws Ø 4.0mm (F.th) x 18mm
217.220	SS	Cannulated Screws Ø 4.0mm (F.th) x 20mm
217.222	SS	Cannulated Screws Ø 4.0mm (F.th) x 22mm
217.224	SS	Cannulated Screws Ø 4.0mm (F.th) x 24mm
217.226	SS	Cannulated Screws Ø 4.0mm (F.th) x 26mm
217.228	SS	Cannulated Screws Ø 4.0mm (F.th) x 28mm
217.230	SS	Cannulated Screws Ø 4.0mm (F.th) x 30mm
217.232	SS	Cannulated Screws Ø 4.0mm (F.th) x 32mm
217.234	SS	Cannulated Screws Ø 4.0mm (F.th) x 34mm
217.235	SS	Cannulated Screws Ø 4.0mm (F.th) x 35mm
217.236	SS	Cannulated Screws Ø 4.0mm (F.th) x 36mm
217.238	SS	Cannulated Screws Ø 4.0mm (F.th) x 38mm
217.240	SS	Cannulated Screws Ø 4.0mm (F.th) x 40mm
217.242	SS	Cannulated Screws Ø 4.0mm (F.th) x 42mm
217.244	SS	Cannulated Screws Ø 4.0mm (F.th) x 44mm
217.245	SS	Cannulated Screws Ø 4.0mm (F.th) x 45mm
217.246	SS	Cannulated Screws Ø 4.0mm (F.th) x 46mm
217.248	SS	Cannulated Screws Ø 4.0mm (F.th) x 48mm
217.250	SS	Cannulated Screws Ø 4.0mm (F.th) x 50mm
217.252	SS	Cannulated Screws Ø 4.0mm (F.th) x 52mm
217.254	SS	Cannulated Screws Ø 4.0mm (F.th) x 54mm
217.255	SS	Cannulated Screws Ø 4.0mm (F.th) x 55mm
217.258	SS	Cannulated Screws Ø 4.0mm (F.th) x 58mm
217.260	SS	Cannulated Screws Ø 4.0mm (F.th) x 60mm
217.265	SS	Cannulated Screws Ø 4.0mm (F.th) x 65mm
217.270	SS	Cannulated Screws Ø 4.0mm (F.th) x 70mm

CANNULATED SCREWS Ø 4.0MM (FULLY THREADED)

217.275	SS	Cannulated Screws Ø 4.0mm (F.th) x 75mm
217.280	SS	Cannulated Screws Ø 4.0mm (F.th) x 80mm
617.230	TT	Cannulated Screws Ø 4.0mm (F.th) x 30mm
617.232	TT	Cannulated Screws Ø 4.0mm (F.th) x 32mm
617.234	TT	Cannulated Screws Ø 4.0mm (F.th) x 34mm
617.235	TT	Cannulated Screws Ø 4.0mm (F.th) x 35mm
617.240	TT	Cannulated Screws Ø 4.0mm (F.th) x 40mm
617.242	TT	Cannulated Screws Ø 4.0mm (F.th) x 42mm
617.245	TT	Cannulated Screws Ø 4.0mm (F.th) x 45mm
617.250	TT	Cannulated Screws Ø 4.0mm (F.th) x 50mm
617.255	TT	Cannulated Screws Ø 4.0mm (F.th) x 55mm
617.260	TT	Cannulated Screws Ø 4.0mm (F.th) x 60mm
617.265	TT	Cannulated Screws Ø 4.0mm (F.th) x 65mm
617.270	TT	Cannulated Screws Ø 4.0mm (F.th) x 70mm
617.275	TT	Cannulated Screws Ø 4.0mm (F.th) x 75mm
617.280	TT	Cannulated Screws Ø 4.0mm (F.th) x 80mm
617.285	TT	Cannulated Screws Ø 4.0mm (F.th) x 85mm
617.290	TT	Cannulated Screws Ø 4.0mm (F.th) x 90mm

Specification

- Available in 10 mm to 90 mm length with 4 mm diameter.
- Non-threaded conical head.
- fully threaded shaft.
- Available in both Titanium and stainless steel.

Indication for use

- Cannulated screws are used as lag screws for fixation of diaphyseal fractures.
- It is used in the fixation of Small Fragment Plates.

