



CAT. NO.	METAL	CANNULATED SCREWS Ø 4.5MM (fully Threaded)
220.716	SS	Cannulated Screws Ø 4.5mm (F.th) x 30mm
220.719	SS	Cannulated Screws Ø 4.5mm (F.th) x 35mm
220.722	SS	Cannulated Screws Ø 4.5mm (F.th) x 40mm
220.725	SS	Cannulated Screws Ø 4.5mm(F.th) x 45mm
220.728	SS	Cannulated Screws Ø 4.5mm (F.th) x 50mm
220.731	SS	Cannulated Screws Ø 4.5mm (F.th) x 55mm
220.734	SS	Cannulated Screws Ø 4.5mm (F.th) x 60mm

Specification

- Available in 30,35,40,45,50,55 & 60mm length with 4.5mm diameter.
- Non-threaded conical head.
- fully threaded shaft.
- Available in stainless steel.

Indication for use

- Cannulated screws are used as lag screws for fixation of diaphyseal fractures.
- It is used in the fixation of Small Fragment Plates.