



CAT. NO.	CARBON HALF RINGS (ADULT)
350.312	CARBON HALF RINGS (ADULT) 80mm ID 120mm OD
350.314	CARBON HALF RINGS (ADULT) 100mm ID 140mm OD
350.315	CARBON HALF RINGS (ADULT) 110mm ID 150mm OD
350.316	CARBON HALF RINGS (ADULT) 120mm ID 160mm OD
350.317	CARBON HALF RINGS (ADULT) 130mm ID 170mm OD
350.318	CARBON HALF RINGS (ADULT) 140mm ID 180mm OD
350.319	CARBON HALF RINGS (ADULT) 150mm ID 190mm OD
350.320	CARBON HALF RINGS (ADULT) 160mm ID 200mm OD
350.322	CARBON HALF RINGS (ADULT) 180mm ID 220mm OD
350.324	CARBON HALF RINGS (ADULT) 200mm ID 240mm OD
350.326	CARBON HALF RINGS (ADULT) 220mm ID 260mm OD
350.328	CARBON HALF RINGS (ADULT) 240mm ID 280mm OD
350.330	CARBON HALF RINGS (ADULT) 260mm ID 300mm OD

Specification

- Carbon half rings (ADULT) available in carbon.
- It is available in 80-200mm I.D & 120-300mm O.D.

Indication for use

- Carbon half rings (ADULT) is available is used for attaching the wires passed through the bones.

