



| CAT. NO. | CARBON HALF RINGS (PEDIATRIC) |
|----------|---|
| 350.412 | CARBON HALF RINGS (PEDIATRIC)80mm ID 112mm OD |
| 350.422 | CARBON HALF RINGS (PEDIATRIC)90mm ID 122mm OD |
| 350.432 | CARBON HALF RINGS (PEDIATRIC) 100mm ID 132mm OD |
| 350.452 | CARBON HALF RINGS (PEDIATRIC)120mm ID 152mm OD |

Specification

- Carbon half rings (PEDIATRIC) available in carbon.
- It is available in 80,90,100,120mm I.D & 112,122,132,152 MM O.D.

Indication for use

- Carbon half rings (PEDIATRIC) is used for attaching the wires passed through the bones.

