



CAT. NO.	METAL	CCS SCREW Ø 7.2mm
703.158	TT	CCS SCREW Ø 7.2mm/24.Th.x 080 mm
703.163	TT	CCS SCREW Ø 7.2mm/24.Th.x 085 mm
703.166	TT	CCS SCREW Ø 7.2mm/24.Th.x 087.5 mm
703.168	TT	CCS SCREW Ø 7.2mm/24.Th.x 090 mm
703.171	TT	CCS SCREW Ø 7.2mm/24.Th.x 092,5 mm
703.173	TT	CCS SCREW Ø 7.2mm/24.Th.x 095 mm

Specification

- Available in 80,85,87,90,92.5 & 95mm length with 7.2mm diameter
- Partially 24mm threaded shaft
- Available in Titanium.

Indication for use

- Cannulated cancellous screws (CCS) are the mainstream internal fixations used for osteosynthesis of femoral neck fractures.