



CAT. NO.	METAL	CORTICAL SCREW 1.2MM STAR HEAD
601.304	TT	CORTICAL SCREW 1.2MM STAR HEAD 4.0MM
601.306	TT	CORTICAL SCREW 1.2MM STAR HEAD 6.0MM
601.308	TT	CORTICAL SCREW 1.2MM STAR HEAD 8.0MM
601.310	TT	CORTICAL SCREW 1.2MM STAR HEAD 10.0MM
601.312	TT	CORTICAL SCREW 1.2MM STAR HEAD 12.0MM
601.314	TT	CORTICAL SCREW 1.2MM STAR HEAD 14.0MM

## Specification

- Self- tapping screw can tap its own hole and hence can ease the process of insertion.
- Compress the plate to the bone and provide interfragmentary compression.
- Available in 4,6,8,10,12 & 14mm length against 1.2mm diameter.
- Non-threaded head.
- Available in Titanium.

## Indication for use

- The cortical screw is indicated for use in hard cortical region of the boneSelf tapping feature eliminates the use of bone tap These screws are also available with Self drilling feature.