



CAT. NO.	METAL	CORTICAL SCREW 2.0MM STAR HEAD
601.356	TT	CORTICAL SCREW 2.0MM STAR HEAD 6.0MM
601.358	TT	CORTICAL SCREW 2.0MM STAR HEAD 8.0MM
601.360	TT	CORTICAL SCREW 2.0MM STAR HEAD 10.0MM
601.362	TT	CORTICAL SCREW 2.0MM STAR HEAD 12.0MM
601.364	TT	CORTICAL SCREW 2.0MM STAR HEAD 14.0MM
601.366	TT	CORTICAL SCREW 2.0MM STAR HEAD 16.0MM
601.368	TT	CORTICAL SCREW 2.0MM STAR HEAD 18.0MM
601.370	TT	CORTICAL SCREW 2.0MM STAR HEAD 20.0MM

Specification

- Self- tapping screw can tap its own hole and hence can ease the process of insertion.
- Compress the plate to the bone and provide interfragmentary compression.
- Available in 6,8,10,12,14,16,18 & 20mm length against 2.0mm diameter.
- Non-threaded head.
- Available in Titanium.

Indication for use

- The cortical screw is indicated for use in hard cortical region of the bone Self tapping feature eliminates the use of bone tap These screws are also available with Self drilling feature.