

CORTICAL SCREWS 2.7 MM SELF TAPPING (STAR HEAD)



CAT. NO. SS	CAT. NO. TT	CORTICAL SCREWS 2.7 MM SELF TAPPING (STAR HEAD)
203.110	603.110	CORTICAL SCREWS 2.7 MM SELF TAPPING 10 MM (STAR HEAD)
203.112	603.112	CORTICAL SCREWS 2.7 MM SELF TAPPING 12 MM (STAR HEAD)
203.114	603.114	CORTICAL SCREWS 2.7 MM SELF TAPPING 14 MM (STAR HEAD)
203.116	603.116	CORTICAL SCREWS 2.7 MM SELF TAPPING 16 MM (STAR HEAD)
203.118	603.118	CORTICAL SCREWS 2.7 MM SELF TAPPING 18 MM (STAR HEAD)
203.120	603.120	CORTICAL SCREWS 2.7 MM SELF TAPPING 20 MM (STAR HEAD)
203.122	603.122	CORTICAL SCREWS 2.7 MM SELF TAPPING 22 MM (STAR HEAD)
203.124	603.124	CORTICAL SCREWS 2.7 MM SELF TAPPING 24 MM (STAR HEAD)
203.126	603.126	CORTICAL SCREWS 2.7 MM SELF TAPPING 26 MM (STAR HEAD)
203.128	603.128	CORTICAL SCREWS 2.7 MM SELF TAPPING 28 MM (STAR HEAD)
203.130	603.130	CORTICAL SCREWS 2.7 MM SELF TAPPING 30 MM (STAR HEAD)
203.132	603.132	CORTICAL SCREWS 2.7 MM SELF TAPPING 32 MM (STAR HEAD)
203.134	603.134	CORTICAL SCREWS 2.7 MM SELF TAPPING 34 MM (STAR HEAD)
203.136	603.136	CORTICAL SCREWS 2.7 MM SELF TAPPING 36 MM (STAR HEAD)
203.138	603.138	CORTICAL SCREWS 2.7 MM SELF TAPPING 38 MM (STAR HEAD)
203.140	603.140	CORTICAL SCREWS 2.7 MM SELF TAPPING 40 MM (STAR HEAD)

Specification

- Self-tapping screw can tap its own hole and hence can ease the process of insertion.
- Compress the plate to the bone and provide inter fragmentary compression.
- Available in 10mm to 40 mm length against 2.7mm diameter.
- Non-threaded head.
- Available in Stainless steel & Titanium.

Indication for use

- The cortical screw is indicated for use in hard cortical region of the bone Self tapping feature eliminates the use of bone tap these screws are also available with self drilling feature.