


CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM



CAT. NO. SS	CAT. NO. TT	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM
206,010	630,010	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 10 MM
206,012	630,012	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 12 MM
206,014	630,014	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 14 MM
206,016	630,016	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 16 MM
206,018	630,018	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 18 MM
206,020	630,020	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 20 MM
206,022	630,022	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 22 MM
206,024	630,024	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 24 MM
206,026	630,026	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 26 MM
206,028	630,028	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 28 MM
206,030	630,030	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 30 MM
206,032	630,032	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 32 MM
206,034	630,034	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 34 MM
206,036	630,036	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 36 MM
206,038	630,038	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 38 MM
206,040	630,040	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 40 MM
206,042	630,042	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 42 MM
206,044	630,044	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 44 MM
206,046	630,046	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 46 MM
206,048	630,048	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 48 MM
206,050	630,050	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 50 MM
206,052	630,052	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 52 MM
206,054	630,054	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 54 MM
206,055	630,055	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 55 MM
206,056	630,056	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 56 MM
206,058	630,058	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 58 MM
206,060	630,060	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 60 MM
206,062	630,062	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 62 MM
206,064	630,064	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 64 MM
206,065	630,065	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 65 MM
206,066	630,066	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 66 MM
206,068	630,068	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 68 MM
206,070	630,070	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 70 MM
206,072	630,072	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 72 MM
206,074	630,074	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 74 MM
206,075	630,075	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 75 MM
206,076	630,076	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 76 MM
206,078	630,078	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 78 MM
206,080	630,080	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 80 MM
206,085	630,085	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 85 MM
206,090	630,090	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 90MM
206,095	630,095	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 95 MM
206,100	630,100	CORTICAL SCREW (SELF TAPPING) Ø 3.5 MM X 100 MM

Specification

- Available in 40mm to 120mm length with 6.5mm diameter
- Partially threaded shaft
- Available in both Stainless Steel & Titanium

Indication for use

- Cancellous screws are indicated for securing plates into metaphyseal and epiphyseal areas of bones
- Deep cuts and widely spaced threads enhances the holding power of the screw in soft bones