


CORTICALSCREWS 2.4 MM SELF TAPPING (STAR HEAD)



CAT. NO. SS	CAT. NO. TT	CORTICALSCREWS 2.4 MM SELF TAPPING (STAR HEAD)
202.050	602.050	CORTICALSCREWS 2.4 MM SELF TAPPING 10 MM (STAR HEAD)
202.052	602.052	CORTICALSCREWS 2.4 MM SELF TAPPING 12 MM (STAR HEAD)
202.054	602.054	CORTICALSCREWS 2.4 MM SELF TAPPING 14 MM (STAR HEAD)
202.056	602.056	CORTICALSCREWS 2.4 MM SELF TAPPING 16 MM (STAR HEAD)
202.058	602.058	CORTICALSCREWS 2.4 MM SELF TAPPING 18 MM (STAR HEAD)
202.060	602.060	CORTICALSCREWS 2.4 MM SELF TAPPING 20 MM (STAR HEAD)
202.062	602.062	CORTICALSCREWS 2.4 MM SELF TAPPING 22 MM (STAR HEAD)
202.064	602.064	CORTICALSCREWS 2.4 MM SELF TAPPING 24 MM (STAR HEAD)
202.066	602.066	CORTICALSCREWS 2.4 MM SELF TAPPING 26 MM (STAR HEAD)
202.068	602.068	CORTICALSCREWS 2.4 MM SELF TAPPING 28 MM (STAR HEAD)
202.070	602.070	CORTICALSCREWS 2.4 MM SELF TAPPING 30 MM (STAR HEAD)
202.072	602.072	CORTICALSCREWS 2.4 MM SELF TAPPING 32 MM (STAR HEAD)
202.074	602.074	CORTICALSCREWS 2.4 MM SELF TAPPING 34 MM (STAR HEAD)
202.076	602.076	CORTICALSCREWS 2.4 MM SELF TAPPING 36 MM (STAR HEAD)
202.078	602.078	CORTICALSCREWS 2.4 MM SELF TAPPING 38 MM (STAR HEAD)
202.080	602.080	CORTICALSCREWS 2.4 MM SELF TAPPING 40 MM (STAR HEAD)

Specification

- Self-tapping screw can tap its own hole and hence can ease the process of insertion.
- It is often fixed at the distal end of the plate.
- Wise-Locking screws favours unicortical screw fixation and thus transfer load to the near cortex.
- Available in 10 mm to 40 mm length against 2.4 mm diameter.
- Threaded Star head.
- Fully threaded shaft.
- Self-tapping tip.
- Available in Stainless steel & Titanium.

Indication for use

- Self-tapping screw can tap its own hole and hence can ease the process of insertion.
- Compress the plate to the bone or create axial compression.