



CAT. NO.	METAL	Cancellous Schanz Screws 6.5MM
233.016	SS	Cancellous Schanz Screws 6.5 mm x 16th x 225mm
233.032	SS	Cancellous Schanz Screws 6.5 mm x 32th x 225mm

Specification

- Cancellous Schanz Screws available in Stainless steel.
- Available in 225mm length with 6.5mm Dia & threaded length 16 & 32.

Indication for use

- Cancellous Schanz Screws to fix the cancellous part of the Tibia & Femur with external fixators.