



CAT. NO. SS	CAT. NO. TT	CANCELLOUS SCREW Ø 4.0 MM (P.TH)
211.010	-	Cancellous Screws Ø 4.0 mm (P.th) x 10 mm
211.012	611.012	Cancellous Screws Ø 4.0 mm (P.th) x 12 mm
211.014	611.014	Cancellous Screws Ø 4.0 mm (P.th) x 14 mm
211.016	611.016	Cancellous Screws Ø 4.0 mm (P.th) x 16 mm
211.018	611.018	Cancellous Screws Ø 4.0 mm (P.th) x 18 mm
211.020	611.020	Cancellous Screws Ø 4.0 mm (P.th) x 20 mm
211.022	-	Cancellous Screws Ø 4.0 mm (P.th) x 22 mm
211.024	611.024	Cancellous Screws Ø 4.0 mm (P.th) x 24 mm
211.025	-	Cancellous Screws Ø 4.0 mm (P.th) x 25 mm
211.026	611.026	Cancellous Screws Ø 4.0 mm (P.th) x 26 mm
211.028	611.028	Cancellous Screws Ø 4.0 mm (P.th) x 28 mm
211.030	611.030	Cancellous Screws Ø 4.0 mm (P.th) x 30 mm
211.032	-	Cancellous Screws Ø 4.0 mm (P.th) x 32 mm
211.034	-	Cancellous Screws Ø 4.0 mm (P.th) x 34 mm
211.035	611.035	Cancellous Screws Ø 4.0 mm (P.th) x 35 mm
211.036	-	Cancellous Screws Ø 4.0 mm (P.th) x 36 mm
211.038	-	Cancellous Screws Ø 4.0 mm (P.th) x 38 mm
211.040	611.040	Cancellous Screws Ø 4.0 mm (P.th) x 40 mm
211.044	-	Cancellous Screws Ø 4.0 mm (P.th) x 44mm
211.045	611.045	Cancellous Screws Ø 4.0 mm (P.th) x 45 mm
211.046	-	Cancellous Screws Ø 4.0 mm (P.th) x 46 mm
211.048	-	Cancellous Screws Ø 4.0 mm (P.th) x 48 mm
211.050	611.050	Cancellous Screws Ø 4.0 mm (P.th) x 50 mm
211.055	-	Cancellous Screws Ø 4.0 mm (P.th) x 55 mm
211.060	-	Cancellous Screws Ø 4.0 mm (P.th) x 60 mm
211.065	611.065	Cancellous Screws Ø 4.0 mm (P.th) x 65 mm
211.070	611.070	Cancellous Screws Ø 4.0 mm (P.th) x 70 mm
211.075	-	Cancellous Screws Ø 4.0 mm (P.th) x 75 mm
211.080	-	Cancellous Screws Ø 4.0 mm (P.th) x 80 mm

## Specification

- Available in 10mm to 80mm length with 4.0mm diameter.
- Partially threaded shaft.
- Available in both Stainless Steel & Titanium

## Indication for use

- A screw with relatively coarser thread and often with a smooth, unthreaded portion, which allows it to act as a lag screw and to anchor in soft medullary bone.

