



CAT. NO.	METAL	Cancellous Screws Ø 4.0 mm (16.th)
212.132	SS	Cancellous Screws Ø 4.0 mm (16.th) x 32 mm
212.134	SS	Cancellous Screws Ø 4.0 mm (16.th) x 34 mm
212.135	SS	Cancellous Screws Ø 4.0 mm (16.th) x 35 mm
212.136	SS	Cancellous Screws Ø 4.0 mm (16.th) x 36 mm
212.138	SS	Cancellous Screws Ø 4.0 mm (16.th) x 38 mm
212.140	SS	Cancellous Screws Ø 4.0 mm (16.th) x 40 mm
212.142	SS	Cancellous Screws Ø 4.0 mm (16.th) x 42 mm
212.144	SS	Cancellous Screws Ø 4.0 mm (16.th) x 44 mm
212.145	SS	Cancellous Screws Ø 4.0 mm (16.th) x 45 mm
212.146	SS	Cancellous Screws Ø 4.0 mm (16.th) x 46 mm
212.148	SS	Cancellous Screws Ø 4.0 mm (16.th) x 48 mm
212.150	SS	Cancellous Screws Ø 4.0 mm (16.th) x 50 mm
212.155	SS	Cancellous Screws Ø 4.0 mm (16.th) x 55 mm
212.160	SS	Cancellous Screws Ø 4.0 mm (16.th) x 60 mm
212.165	SS	Cancellous Screws Ø 4.0 mm (16.th) x 65 mm
212.170	SS	Cancellous Screws Ø 4.0 mm (16.th) x 70 mm
212.175	SS	Cancellous Screws Ø 4.0 mm (16.th) x 75 mm

Specification

- Available in 32-75mm length with 4.0mm diameter
- Partially 16mm threaded shaft
- Available in Stainless Steel.

Indication for use

- Cancellous screws are indicated for securing plates into metaphyseal and epiphyseal areas of bones.
- Deep cuts and widely spaced threads enhance the holding power of the screw in soft bones.