



CAT. NO.	METAL	Cannulated Screws Ø 3.5mm (F.th)
217.810	SS	Cannulated Screws Ø 3.5mm (F.th) x 10mm
217.812	SS	Cannulated Screws Ø 3.5mm (F.th) x 12mm
217.814	SS	Cannulated Screws Ø 3.5mm (F.th) x 14mm
217.816	SS	Cannulated Screws Ø 3.5mm (F.th) x 16mm
217.818	SS	Cannulated Screws Ø 3.5mm (F.th) x 18mm
217.820	SS	Cannulated Screws Ø 3.5mm (F.th) x 20mm

## Specification

- Available in 10,12,14,16,18 & 20mm length with 3.5mm diameter.
- Non-threaded conical head.
- fully threaded shaft.
- Available in stainless steel.

## Indication for use

- 3.5mm (F.th) Cannulated Screws are intended for the fixation of fractures, osteotomies, fusions and Nonunions of small bones and small bone fragments.

