



CAT. NO. SS	Condylar Blade Plate 95° for Adults
302.150	Condylar Plate 95° - 05 Holes, Blade Length- 50 mm
302.160	Condylar Plate 95° - 05 Holes, Blade Length- 60 mm
302.170	Condylar Plate 95° - 05 Holes, Blade Length- 70 mm
302.180	Condylar Plate 95° - 05 Holes, Blade Length- 80 mm
302.350	Condylar Plate 95° - 07 Holes, Blade Length- 50 mm
302.360	Condylar Plate 95° - 07 Holes, Blade Length- 60 mm
302.370	Condylar Plate 95° - 07 Holes, Blade Length- 70 mm
302.380	Condylar Plate 95° - 07 Holes, Blade Length- 80 mm
302.450	Condylar Plate 95° - 09 Holes, Blade Length- 50 mm
302.460	Condylar Plate 95° - 09 Holes, Blade Length- 60 mm
302.470	Condylar Plate 95° - 09 Holes, Blade Length- 70 mm
302.480	Condylar Plate 95° - 09 Holes, Blade Length- 80 mm
302.650	Condylar Plate 95° - 12 Holes, Blade Length- 50 mm
302.660	Condylar Plate 95° - 12 Holes, Blade Length- 60 mm
302.670	Condylar Plate 95° - 12 Holes, Blade Length- 70 mm
302.680	Condylar Plate 95° - 12 Holes, Blade Length- 80 mm
302.850	Condylar Plate 95° - 14 Holes, Blade Length- 50 mm
302.860	Condylar Plate 95° - 14 Holes, Blade Length- 60 mm
302.870	Condylar Plate 95° - 14 Holes, Blade Length- 70 mm
302.880	Condylar Plate 95° - 14 Holes, Blade Length- 80 mm
302.890	Condylar Plate 95° - 16 Holes, Blade Length- 50 mm
302.891	Condylar Plate 95° - 16 Holes, Blade Length- 60 mm
302.892	Condylar Plate 95° - 16 Holes, Blade Length- 70 mm
302.893	Condylar Plate 95° - 16 Holes, Blade Length- 80 mm
302.950	Condylar Plate 95° - 16 Holes, Blade Length- 50 mm
302.960	Condylar Plate 95° - 16 Holes, Blade Length- 60 mm
302.970	Condylar Plate 95° - 16 Holes, Blade Length- 70 mm
302.980	Condylar Plate 95° - 16 Holes, Blade Length- 80 mm
302.995	Condylar Plate 95° - 18 Holes, Blade Length- 50 mm
302.996	Condylar Plate 95° - 18 Holes, Blade Length- 60 mm
302.997	Condylar Plate 95° - 18 Holes, Blade Length- 70 mm
302.998	Condylar Plate 95° - 18 Holes, Blade Length- 80 mm

Specification

- The angle between blade and shaft is 95°.
- It uses 4.5/6.5mm Cortical screws Cancellous screws.
- Two holes next to the blade are suitable for 6.5mm Cancellous bone screw.
- Plates available with Blade 50,60,70 & 80mm size.

Indication for use

- Condylar Blade Plate As buttress for the femoral condyles with multiple fragment fractures.