

CONDYLAR BUTTRESS PLATE 4.5 FOR LEFT FEMORAL CONDYLES



CAT. NO. SS	Condylar Buttress Plate 4.5 for Left Femoral Condyles
303.005L	Condylar Buttress Plate - 05 Holes- Left
303.005R	Condylar Buttress Plate - 05 Holes- Right
303.006L	Condylar Buttress Plate - 06 Holes- Left
303.006R	Condylar Buttress Plate - 06 Holes- Right
303.007L	Condylar Buttress Plate - 07 Holes- Left
303.007R	Condylar Buttress Plate - 07 Holes- Right
303.008L	Condylar Buttress Plate - 08 Holes- Left
303.008R	Condylar Buttress Plate - 08 Holes- Right
303.009L	Condylar Buttress Plate - 09 Holes- Left
303.009R	Condylar Buttress Plate - 09 Holes- Right
303.010L	Condylar Buttress Plate - 10 Holes- Left
303.010R	Condylar Buttress Plate - 10 Holes- Right
303.011L	Condylar Buttress Plate - 11 Holes- Left
303.011R	Condylar Buttress Plate - 11 Holes- Right
303.012L	Condylar Buttress Plate - 12 Holes- Left
303.012R	Condylar Buttress Plate - 12 Holes- Right
303.014L	Condylar Buttress Plate - 14 Holes- Left
303.014R	Condylar Buttress Plate - 14 Holes- Right

Specification

- Plates available in Stainless steel.
- Plates available with 5,6,7,8,9,10,11,12 & 14 holes.
- It accepts 4.5/6.5mm Cortical screws Cancellous screws.

Indication for use

- Condylar Buttress Plate use in six fresh-frozen human knee specimens, segmental osteotomies were created to mimic supracondylar femoral fractures. The medial cortex was completely removed to make the fracture unstable to varus deformity. The fracture was fixed with a lateral condylar buttress plate using 4.5 mm screws.