



| CAT. NO. | METAL | Connecting Rod                  |
|----------|-------|---------------------------------|
| 368.250  | SS    | Connecting Rod Ø 2.5mm x 75 mm  |
| 368.252  | SS    | Connecting Rod Ø 2.5mm x 100 mm |
| 368.254  | SS    | Connecting Rod Ø 2.5mm x 115 mm |
| 368.256  | SS    | Connecting Rod Ø 2.5mm x 125 mm |
| 368.258  | SS    | Connecting Rod Ø 2.5mm x 150 mm |
| 368.260  | SS    | Connecting Rod Ø 2.5mm x 175 mm |
| 368.262  | SS    | Connecting Rod Ø 2.5mm x 200 mm |
| 368.264  | SS    | Connecting Rod Ø 2.5mm x 225 mm |
| 368.266  | SS    | Connecting Rod Ø 2.5mm x 250 mm |
| 368.268  | SS    | Connecting Rod Ø 2.5mm x 275 mm |
| 368.270  | SS    | Connecting Rod Ø 2.5mm x 300 mm |
| 368.325  | SS    | Connecting Rod Ø 3.0mm x 75 mm  |
| 368.327  | SS    | Connecting Rod Ø 3.0mm x 100 mm |
| 368.329  | SS    | Connecting Rod Ø 3.0mm x 115 mm |
| 368.331  | SS    | Connecting Rod Ø 3.0mm x 125 mm |
| 368.333  | SS    | Connecting Rod Ø 3.0mm x 150 mm |
| 368.335  | SS    | Connecting Rod Ø 3.0mm x 175 mm |
| 368.337  | SS    | Connecting Rod Ø 3.0mm x 200 mm |
| 368.339  | SS    | Connecting Rod Ø 3.0mm x 225 mm |
| 368.341  | SS    | Connecting Rod Ø 3.0mm x 250 mm |
| 368.343  | SS    | Connecting Rod Ø 3.0mm x 275 mm |
| 368.345  | SS    | Connecting Rod Ø 3.0mm x 300 mm |
| 368.352  | SS    | Connecting Rod Ø 3.5mm x 75 mm  |
| 368.354  | SS    | Connecting Rod Ø 3.5mm x 100 mm |
| 368.356  | SS    | Connecting Rod Ø 3.5mm x 115 mm |
| 368.358  | SS    | Connecting Rod Ø 3.5mm x 125 mm |
| 368.360  | SS    | Connecting Rod Ø 3.5mm x 150 mm |
| 368.362  | SS    | Connecting Rod Ø 3.5mm x 175 mm |
| 368.364  | SS    | Connecting Rod Ø 3.5mm x 200 mm |
| 368.366  | SS    | Connecting Rod Ø 3.5mm x 225 mm |
| 368.368  | SS    | Connecting Rod Ø 3.5mm x 250 mm |
| 368.370  | SS    | Connecting Rod Ø 3.5mm x 275 mm |
| 368.372  | SS    | Connecting Rod Ø 3.5mm x 300 mm |

## Specification

- Connecting Rod available in Stainless steel.
- Available length 75-300mm with 2.5, 3.0 & 3.5mm dia.

## Indication for use

- Orthopaedic Connecting Rods are widely utilized for stabilizing the soft and bone tissues. These are simple to fix and boast of strong construction. The provided rods are of precise dimension. The rods allow for external fixation and are suited for many surgical treatments.

