



CAT. NO. SS	CAT. NO. TT	DCS BARREL PLATE 95° WITH COMBI HOLES
319.606	619.606	DCS BARREL PLATE 95° WITH COMBI HOLES - 06 HOLES, L - 100 MM
319.608	619.608	DCS BARREL PLATE 95° WITH COMBI HOLES - 08 HOLES, L - 130 MM
319.610	619.610	DCS BARREL PLATE 95° WITH COMBI HOLES - 10 HOLES, L - 163 MM
319.612	619.612	DCS BARREL PLATE 95° WITH COMBI HOLES - 12 HOLES, L - 198 MM

Specification

- Plates available with 6,8,10 & 12 holes.
- Shaft of the plate mate with 4.5mm Cortical screw, Lcp 5.0mm and barrel mate with 12.5mm DHS/DCS Screws (with Compression Screws)

Indication for use

It is indicated to fix following fractures of the distal femur -

- Intercondylar fractures
- Supracondylar fractures
- Unicondylar fractures

