

## DHS BARREL PLATE WITH DCP HOLE-AO TYPE, 125°- LONG (BARREL LENGTH- 38MM)



CAT. NO. SS	DHS Barrel Plate with DCP Hole -AO type, 125° - Long (Barrel length - 38mm)
289.003	DHS L/B Plate with DCP Hole 125° -03 HOLES, L - 62 MM
289.004	DHS L/B Plate with DCP Hole 125° -04 HOLES, L - 78 MM
289.005	DHS L/B Plate with DCP Hole 125° -05 HOLES, L - 94 MM
289.006	DHS L/B Plate with DCP Hole 125° -06 HOLES, L - 110 MM

### Specification

- Plates available in Stainless steel.
- Plates available with 3,4,5 & 6 holes.
- Plate accepts 4.5mm Cortical screw and barrel mate with 12.5mm DHS/DCS Screws (with Compression Screws).

### Indication for use

- The Dynamic Hip Plate offered by us has been designed for providing strong and stable internal fixation for various applications in intertrochanteric, subtrochanteric and basilar neck fractures, while causing minimal soft tissue irritation.