

DHS BARREL PLATE WITH DCP HOLE AO TYPE, 130°- LONG (BARREL LENGTH- 38MM)



CAT. NO. SS	CAT. NO. TT	DHS Barrel Plate with DCP Hole AO type, 130° - Long (Barrel length - 38mm)
289.102	689.102	DHS L/B Plate with DCP Hole 130° - 02 Holes, L- 46mm
289.103	689.103	DHS L/B Plate with DCP Hole 130° - 03 Holes, L- 62mm
289.104	689.104	DHS L/B Plate with DCP Hole 130° - 04 Holes, L- 78mm
289.105	689.105	DHS L/B Plate with DCP Hole 130° - 05 Holes, L- 94mm
289.106	689.106	DHS L/B Plate with DCP Hole 130° - 06 Holes, L- 110mm
289.107	689.107	DHS L/B Plate with DCP Hole 130° - 07 Holes, L- 126mm
289.108	689.108	DHS L/B Plate with DCP Hole 130° - 08 Holes, L- 142mm
289.110	689.110	DHS L/B Plate with DCP Hole 130° - 10 Holes, L- 174mm

Specification

- Plates available in Stainless steel & Titanium.
- Plates available with 2,3,4,5,6,7,8 & 10 holes.
- Plate accepts 4.5mm Cortical screw and barrel mate with 12.5mm DHS/DCS Screws (with Compression Screws).

Indication for use

- The Dynamic Hip Plate offered by us has been designed for providing strong and stable internal fixation for various applications in intertrochanteric, subtrochanteric and basilar neck fractures, while causing minimal soft tissue irritation.