

DHS BARREL PLATE WITH DCP HOLE AO TYPE, 135°- LONG (BARREL LENGTH- 38MM)



CAT. NO. SS	CAT. NO. TT	DHS Barrel Plate with DCP Hole AO type, 13 5°- Long (Barrel length - 38mm)
289.202	689.202	DHS L/B Plate with DCP Hole 135° - 02 Holes, L- 46mm
289.203	689.203	DHS L/B Plate with DCP Hole 135° - 03 Holes, L- 62mm
289.204	689.204	DHS L/B Plate with DCP Hole 135° - 04 Holes, L- 78mm
289.205	689.205	DHS L/B Plate with DCP Hole 135° - 05 Holes, L- 94mm
289.206	-	DHS L/B Plate with DCP Hole 135° - 06 Holes, L- 110mm
289.207	-	DHS L/B Plate with DCP Hole 135° - 07 Holes, L- 126mm
289.208	689.208	DHS L/B Plate with DCP Hole 135° - 08 Holes, L- 142mm
289.209	-	DHS L/B Plate with DCP Hole 135° - 09 Holes, L- 142mm
289.210	689.210	DHS L/B Plate with DCP Hole 135° - 10 Holes, L- 174mm

Specification

- Plates available in Stainless steel & Titanium.
- Plates available with 2,3,4,5,6,7,8 & 10 holes.
- Plate accepts 4.5mm Cortical screw and barrel mate with 12.5mm DHS/DCS Screws (with Compression Screws).

Indication for use

- The Dynamic Hip Plate offered by us has been designed for providing strong and stable internal fixation for various applications in intertrochanteric, subtrochanteric and basilar neck fractures, while causing minimal soft tissue irritation.