

DHS BARREL PLATE WITH DCP HOLE AO TYPE, 150°- LONG (BARREL LENGTH- 38MM)



CAT. NO.	DHS Barrel Plate with DCP Hole AO type, 1 50° - Long (Barrel length - 38mm)
SS	
289.402	DHS L/B Plate with DCP Hole 1 50° - 02 Holes, L - 46mm
289.404	DHS L/B Plate with DCP Hole 1 50° - 04 Holes, L - 78mm
289.405	DHS L/B Plate with DCP Hole 1 50° - 05 Holes, L - 94mm
289.406	DHS L/B Plate with DCP Hole 1 50° - 06 Holes, L - 110mm

Specification

- Plates available in Stainless steel.
- Plates available with 2,4,5 & 6 holes.
- Plate accepts 4.5mm Cortical screw and barrel mate with 12.5mm DHS/DCS Screws (with Compression Screws).

Indication for use

- The Dynamic Hip Plate offered by us has been designed for providing strong and stable internal fixation for various applications in intertrochanteric, subtrochanteric and basilar neck fractures, while causing minimal soft tissue irritation.