



CAT. NO. SS	CAT. NO. TT	DHS L/B PLATE WITH COMBI HOLE 135°
319.502	619.502	DHS L/B PLATE WITH COMBI HOLE 135° - 02 HOLES
319.503	619.503	DHS L/B PLATE WITH COMBI HOLE 135° - 03 HOLES
319.504	619.504	DHS L/B PLATE WITH COMBI HOLE 135° - 04 HOLES, L - 78 MM
319.505	619.505	DHS L/B PLATE WITH COMBI HOLE 135° - 05 HOLES, L - 94 MM
319.506	619.506	DHS L/B PLATE WITH COMBI HOLE 135° - 06 HOLES, L - 110MM
319.507	619.507	DHS L/B PLATE WITH COMBI HOLE 135° - 07 HOLES, L - 126MM
319.508	619.508	DHS L/B PLATE WITH COMBI HOLE 135° - 08 HOLES, L - 142MM
319.509	619.509	DHS L/B PLATE WITH COMBI HOLE 135° - 09 HOLES, L - 158MM
319.510	619.510	DHS L/B PLATE WITH COMBI HOLE 135° - 10 HOLES, L - 174MM

Specification

- Available in both Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5mm / 5.0mm cortical and locking screws.
- Plate available with 2,3,4,5,6,7,8,9 &10 holes.

Indication for use

- Wise Lock DHS Plate is intended for use in fixation of fractures to the proximal femur.
- The system is indicated for use in trochanteric, pertrochanteric, intertrochanteric, and basilar neck fracture.