

## DHS L/B PLATE WITH COMBI HOLE 130° (A.O TYPE, LONG)



CAT. NO. SS	CAT. NO. TT	DHS L/B Plate with Combi Hole 130° (A.O Type, Long)
-	619.532	DHS L/B Plate with Combi Hole 130° - 02 Holes
319.533	619.533	DHS L/B Plate with Combi Hole 130° - 03 Holes
319.534	619.534	DHS L/B Plate with Combi Hole 130° - 04 Holes
319.535	619.535	DHS L/B Plate with Combi Hole 130° - 05 Holes
319.536	619.536	DHS L/B Plate with Combi Hole 130° - 06 Holes
319.538	-	DHS L/B Plate with Combi Hole 130° - 08 Holes
319.540	-	DHS L/B Plate with Combi Hole 130° - 10 Holes

### Specification

- Available in both Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5mm/5.0mm cortical and locking screws & DHS screws.
- Plate available with 2,3,4,5,6,7,8,9 & 10 holes.

### Indication for use

#### It is indicated to fix-

- fractures of the proximal femur.
- trochanteric, pertrochanteric.
- Basilar neck fracture.

