

## DHS L/B PLATE WITH COMBI HOLE 125° (A.O TYPE, LONG)



CAT. NO.	DHS L/B Plate with Combi Hole 125° (A.O Type, Long)
TT	
619.564	DHS L/B Plate with Combi Hole 125° - 4 Holes

### Specification

- Available in Titanium.
- Plate features locking compression holes and round locking holes that accepts 4.5mm/5.0mm cortical and locking screws.
- Plate available with 4 holes.

### Indication for use

- It is indicated to fix-
  - fractures of the proximal femur.
  - trochanteric, pertrochanteric.
  - Basilar neck fracture.