



CAT. NO. SS	DHS L/B Plate with Combi Hole 125°
319,562	DHS L/B Plate with Combi Hole 125° - 02 Holes
319,563	DHS L/B Plate with Combi Hole 125° - 03 Holes
319,564	DHS L/B Plate with Combi Hole 125° - 04 Holes
319,565	DHS L/B Plate with Combi Hole 125° - 05 Holes
319,566	DHS L/B Plate with Combi Hole 125° - 06 Holes
319,568	DHS L/B Plate with Combi Hole 125° - 08 Holes
319,569	DHS L/B Plate with Combi Hole 125° - 09 Holes
319,570	DHS L/B Plate with Combi Hole 125° - 10 Holes

## Specification

- Available in Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5mm/5.0mm cortical and locking screws.
- Plate available with 2,3,4,5,6,8,9 & 10 holes.

## Indication for use

### It is indicated to fix-

- fractures of the proximal femur.
- trochanteric, pertrochanteric.
- Basilar neck fracture.