

DHS L/B PLATE WITH COMBI HOLE 140° (A.O TYPE, LONG)



CAT. NO. TT	DHS L/B Plate with Combi Hole 140° (A.O Type, Long)
619.804	DHS L/B Plate with Combi Hole 1 40° - 4 Holes

Specification

- Available in Titanium.
- Plate features locking compression holes and round locking holes that accepts 4.5mm/5.0mm cortical and locking screws.
- Plate available with 4 holes.

Indication for use

- It is indicated to fix-
 - fractures of the proximal femur.
 - trochanteric, pertrochanteric.
 - Basilar neck fracture.

