



DHS L/B Plate with DCP Hole 145°	
	DHS L/B Plate with DCP Hole 145° - 04 Holes, L- 78 mm
	DHS L/B Plate with DCP Hole 145° - 05 Holes, L- 94 mm
	DHS L/B Plate with DCP Hole 145° - 06 Holes, L- 94 mm
	DHS L/B Plate with DCP Hole 145° - 03 Holes, L- 62 mm
	DHS L/B Plate with DCP Hole 145° - 04 Holes, L- 78 mm
	DHS L/B Plate with DCP Hole 145° - 05 Holes, L- 94 mm
	DHS SLB Plate with DCP Hole 145° - 06 Holes, L- 110 mm

Specification

- Available in Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5mm/5.0mm cortical and locking screws.
- Plate available with 3,4,5 & 6 holes.

Indication for use

It is indicated to fix -

- fractures of the proximal femur
- trochanteric, pertrochanteric
- Basilar neck fracture

