



CAT. NO. SS	CAT. NO. TT	DHS S/B Plate with Combi Hole 130°
319.723	619.723	DHS S/B Plate with Combi Hole 130° - 03 Holes
319.724	619.724	DHS S/B Plate with Combi Hole 130° - 04 Holes, L - 78 mm
319.725	619.725	DHS S/B Plate with Combi Hole 130° - 05 Holes, L - 94 mm
319.726	619.726	DHS S/B Plate with Combi Hole 130° - 06 Holes, L - 110 mm
319.728	619.728	DHS S/B Plate with Combi Hole 130° - 08 Holes, L - 142 mm
319.730	619.730	DHS S/B Plate with Combi Hole 130° - 10 Holes, L - 174 mm

Specification

- Available in both Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5mm/5.0mm cortical and locking screws & DHS screws.
- Plate available with 3,4,5,6,8 & 10 holes.

Indication for use

It is indicated to fix-

- fractures of the proximal femur
- trochanteric, pertrochanteric, intertrochanteric, and basilar neck fracture