



CAT. NO. SS	CAT. NO. TT	DHS S/B Plate with Combi Hole 135°
319.703	619.703	DHS S/B Plate with Combi Hole 135° - 03 Holes,
319.704	619.704	DHS S/B Plate with Combi Hole 135° - 04 Holes, L- 78 mm
319.705	619.705	DHS S/B Plate with Combi Hole 135° - 05 Holes, L- 94 mm
319.706	619.706	DHS S/B Plate with Combi Hole 135° - 06 Holes, L- 110mm
319.708	619.708	DHS S/B Plate with Combi Hole 135° - 08 Holes, L- 142mm
319.710	619.710	DHS S/B Plate with Combi Hole 135° - 10 Holes, L- 174mm

Specification

- Available in both Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5mm/5.0mm cortical and locking screws & DHS screws.
- Plate available with 3,4,5,6,8 & 10 holes.

Indication for use

- It is indicated to fix-
 - fractures of the proximal femur
 - trochanteric, pertrochanteric, intertrochanteric, and basilar neck fracture

