

DHS SCREWS WITH COMPRESSION SCREW (STANDARD THREAD)



CAT. NO. SS	CAT. NO. TT	DHS Screws with Compression Screw (Standard Thread)
221.001	-	Spare Compression Screw for DHS Screws
221.050	621.050	DHS Screws with Compression Screw, L - 050 mm
221.055	621.055	DHS Screws with Compression Screw, L - 055 mm
221.060	621.060	DHS Screws with Compression Screw, L - 060 mm
221.065	621.065	DHS Screws with Compression Screw, L - 065 mm
221.070	621.070	DHS Screws with Compression Screw, L - 070 mm
221.075	621.075	DHS Screws with Compression Screw, L - 075 mm
221.080	621.080	DHS Screws with Compression Screw, L - 080 mm
221.085	621.085	DHS Screws with Compression Screw, L - 085 mm
221.090	621.090	DHS Screws with Compression Screw, L - 090 mm
221.095	621.095	DHS Screws with Compression Screw, L - 095 mm
221.100	621.100	DHS Screws with Compression Screw, L - 100 mm
221.105	621.105	DHS Screws with Compression Screw, L - 105 mm
221.110	621.110	DHS Screws with Compression Screw, L - 110 mm
221.115	-	DHS Screws with Compression Screw, L - 115 mm

Specification

- Available in 50 mm to 115 mm length
- Standard threaded shaft.
- Available in both Titanium and stainless steel.

Indication for use

- DHS screw Use with DHS plate trochanteric & fracture neck femur. Take less space in fracture neck femur compare to standard DHS screw thus provides more space for additional fixation to provide two Point fixation.

