



CAT. NO.	METAL	Danham Pin (Centrally Threaded)
240.245	SS	Danham Pin Ø 4.5 mm x 225mm
240.245S	SS	Danham Pin Ø 4.5 mm x 225mm (Sterilized)
240.254	SS	Danham Pin Ø 5.0 mm x 175mm
240.255	SS	Danham Pin Ø 5.0 mm x 225mm
240.256	SS	DanhamPin Ø 6.0 mm x 175mm

## Specification

- Danham Pin (Centrally Threaded) available in Stainless steel.
- Available length 175 & 225mm with 4.5mm dia.

## Indication for use

- Danham Pin used for Fracture Fixation.