

DYNAMIC COMPRESSION PLATES (LC-DCP NARROW) FOR 4.5 MM SCREWS



CAT. NO. SS	CAT. NO. TT	Dynamic Compression Plates (LC -DCP Narrow) for 4.5 mm Screws
285.206	778.106	LC-DCP Narrow for 4.5 mm Screws - 06 Holes, L - 106 mm
285.207	778.107	LC-DCP Narrow for 4.5 mm Screws - 07 Holes, L - 124 mm
285.208	778.108	LC-DCP Narrow for 4.5 mm Screws - 08 Holes, L - 142 mm
285.209	778.109	LC-DCP Narrow for 4.5 mm Screws - 09 Holes, L - 160 mm
285.210	778.110	LC-DCP Narrow for 4.5 mm Screws - 10 Holes, L - 178 mm
285.211		LC-DCP Narrow for 4.5 mm Screws - 11 Holes, L - 196 mm
285.212		LC-DCP Narrow for 4.5 mm Screws - 12 Holes, L - 214 mm
285.214		LC-DCP Narrow for 4.5 mm Screws - 14 Holes, L - 250 mm
285.216		LC-DCP Narrow for 4.5 mm Screws - 16 Holes, L - 250 mm

Specification

- Plates available in Stainless steel.
- Plates available with 3,4,5 & 6 holes.
- Plate accepts 4.5mm Cortical screw and barrel mate with 12.5mm DHS/DCS Screws (with Compression Screws).

Indication for use

- Narrow Dynamic Compression Plates are indicated for fixation of various long bones, such as the humerus, femur and tibia.
- These are also used in fixation of peri-prosthetic fractures, osteopenic bone and fixation of nonunions and malunions in adult patients. The oval shaped hole is designed to create dynamic pressure between the fractured bone fragments.